

A top-down view of a hydrotherapy tub, illuminated with a vibrant blue light. The tub's interior is visible, showing a central human silhouette and several rows of small, glowing blue lights. The tub is set against a dark background, with a bright blue vertical bar on the right side.

HYDROTHERAPY





sales@btlnet.com
www.btlnet.com

All rights reserved. Although every care has been taken to provide accurate and up-to-date information, no responsibility can be accepted for any mistakes or misprints contained within this catalogue. The products and their specifications, as shown in this catalogue, may change without notice.



HYDROTHERAPY

The new series of BTL whirlpools combines clean, elegant shapes with high-quality materials and perfect workmanship. The word “simplicity” best characterizes their design and operation. The unique features, such as the colour touch screen user-interface or the one-touch start of therapy programmes make operation extremely fast and easy. A soft blue glow gently highlights the exterior base of the baths. The baths are designed for all kinds of professional facilities. The series consists of luxurious baths for wellness centres and total-immersion whirlpools in two sizes.

A white ceramic bowl is shown in a close-up, low-angle shot. The bowl is illuminated from below, creating a bright blue glow that highlights its rim and the surface it sits on. The background is dark, making the white bowl and the blue light stand out. The text "AWESOME LOOK AWESOME" is overlaid on the left side of the bowl.


AWESOME LOOK
AWESOME



PERFORMANCE

OMEGA

UNBEATABLE COMFORT & LUXURY



Millions of
bubbles for
ultimate
relaxation



Music option:
Built-in sound
system



Up to 266
water jets



Aromatherapy and
bath additives



Colour therapy



Pre-set protocols
with 'Hot-keys'

RELAXATION BATHS FOR SPAS & WELLNESS CENTRES

SUPERB ERGONOMICS WITH OPTIMALLY-POSITIONED MASSAGE JETS



MAIN FEATURES & OPTIONS

- ◇ Comfortable operation using the colour touch display
- ◇ Ergonomic underwater massage
- ◇ Multi-zone therapy system
- ◇ 'Hot-keys' for frequently used therapies
- ◇ Pre-set therapy protocols
- ◇ Suitable for different kinds of water: Fresh, salt, carbonated or re-circulated
- ◇ Automatic disinfection programme
- ◇ High-quality acrylate bath surface for easy cleaning
- ◇ LED lighting of bath exterior
- ◇ Developed for reliability and long life

HYDROMASSAGE

- ◇ Up to 266 water jets
- ◇ Up to 8 hydromassage zones:
 - Soles of the feet
 - Lower limbs
 - Outer calves and knees
 - Inner calves and knees
 - Buttocks
 - Outer thighs and hips
 - Lower and upper back
 - Upper limbs and shoulders

Model	Omega 10 Deluxe	Omega 20 Deluxe	Omega 30 Deluxe
Specifications			
HYDROMASSAGE			
Number of hydromassage jets	104	168	266
Hydromassage of the following body zones:	2 zones	4 zones	8 zones
Soles of the feet			1
Lower limbs		1	2
Outer calves and knees		2	3
Inner calves and knees	1		4
Buttocks		3	5
Outer thighs and hips			6
Lower and upper back	2		7
Upper limbs and shoulders		4	8
Regulation of hydromassage intensity	●	●	●
Hydromassage with added air	●	●	●
UNDERWATER MASSAGE			
	○	○	○
STANDARD SPECIFICATIONS & OPTIONS			
LCD touch display 5.7"	●	●	●
Water filling temperature indicator	●	●	●
Bath temperature indicator	●	●	●
Automatic bath filling and emptying	●	●	●
Water level indicator with acoustic signaling	●	●	●
Thermostatic mixing tap	●	●	●
Slide-in shower	●	●	●
'Hot-keys' for frequently used programmes	●	●	●
Colour therapy option (Chromatherapy)	○	○	●
Music option (Built-in sound system)	○	○	○
Automatic disinfection programme	○	○	○
Carbonated water bath (CO ₂ saturator)	○	○	○
Saltwater option / re-circulated water option	○	○	○
Foot support	○	○	●
External lighting	○	○	●
TECHNICAL SPECIFICATIONS			
Utility bath volume	220 l		
Dimensions: Length / Width / Height	2370 / 890 / 770 mm		
Main power supply	230 V / 50 Hz		

● – standard ○ – optional



Omega 10
Deluxe



Omega 20
Deluxe



Omega 30
Deluxe



MULTI-ZONE THERAPY SYSTEM FOR COMFORTABLE MASSAGE



HYDROTHERAPY



ONE-TOUCH STA

RT




ALFA

FULL-IMMERSION BATHS FOR ALL TYPES OF MEDICAL FACILITIES & SPAS




Hydromassage
jets



Air massage jets



Ergonomic
underwater
massage



Carbonated
water massage



Aromatherapy
and bath
additives



Colour therapy

FOR ALL THERAPEUTIC AND RESTORATIVE PROCEDURES

POWERFUL MASSAGE SYSTEMS



MAIN FEATURES & OPTIONS

- ◇ Comfortable operation using the colour touch display or buttons
- ◇ Ergonomic underwater massage
- ◇ Multi-zone therapy system
- ◇ 'Hot-keys' for frequently used therapies
- ◇ Pre-set therapy protocols
- ◇ Suitable for different kinds of water: Fresh, salt, carbonated or re-circulated
- ◇ Automatic disinfection programme
- ◇ High-quality acrylate bath surface for easy cleaning
- ◇ LED lighting of bath exterior
- ◇ Developed for reliability and long life

HYDROMASSAGE


- ◇ Up to 8 water jets and 30 micro-jets
- ◇ Up to 4 hydromassage zones:
 - Soles of the feet
 - Lower limbs, outer calves and knees
 - Buttocks, thighs and hips
 - Back, upper limbs and shoulders

AIR MASSAGE

- ◇ Up to 122 air massage jets
- ◇ Up to 2 air massage zones:
 - Lower limbs
 - Hips and buttocks

Model	ALFA 10	ALFA 20	ALFA 30	ALFA 40	ALFA 50	ALFA 60	ALFA 70
Specifications			ALFA 30 Deluxe	ALFA 40 Deluxe		ALFA 60 Deluxe	ALFA 70 Deluxe
HYDROMASSAGE							
Number of hydromassage jets			4	8		8	8
Number of hydromassage micro-jets			12	24		24	30
Hydromassage of the following body zones:			4 zones	4 zones		4 zones	4 zones
Soles of the feet			1	1		1	1
Lower limbs			2	2		2	2
Outer calves and knees			3	3		3	3
Buttocks			3	3		3	3
Outer thighs and hips			4	4		4	4
Lower and upper back			4	4		4	4
Upper limbs and shoulders							
Regulation of hydromassage intensity			●	●		●	●
Hydromassage with added air			●	●		●	●
AIR MASSAGE							
Number of air massage jets					122	104	122
Air massage of the following body zones:					1 zone	2 zones	2 zones
Lower limbs					1	1	1
Hips and buttocks						2	2
Regulation of air massage intensity					●	●	●
UNDERWATER MASSAGE							
	○	●	○	○	○	○	○
CARBONATED WATER BATH (CO₂ saturator)							
	●	○	○	○	○	○	○
STANDARD SPECIFICATIONS & OPTIONS							
Water level indicator with acoustic signaling	●	●	●	●	●	●	●
High water temperature indicator	●	●	●	●	●	●	●
Slide-in shower	●	●	●	●	●	●	●
Thermostatic mixing tap	○	○	○	○	○	●	●
Saltwater option /re-circulated water option	○	○	○	○	○	○	○
Foot support	○	○	○	○	○	○	○
External lighting	○	○	○	○	○	○	○
DELUXE LINE SPECIFICATIONS & OPTIONS							
LCD touch display 5.7"			● (Deluxe)	● (Deluxe)		●	●
Water filling & bath temperature indicator			● (Deluxe)	● (Deluxe)		●	●
Automatic bath filling and emptying			● (Deluxe)	● (Deluxe)		●	●
Thermostatic mixing tap			● (Deluxe)	● (Deluxe)		●	●
'Hot-keys' for frequently used programmes			● (Deluxe)	● (Deluxe)		●	●
Automatic drying of air massage jets						●	●
Colour therapy option (Chromatherapy)			○ (Deluxe)	○ (Deluxe)		○	○
Music option (Built-in sound system)			○ (Deluxe)	○ (Deluxe)		○	○
Automatic disinfection programme			○ (Deluxe)	○ (Deluxe)		○	○
TECHNICAL SPECIFICATIONS							
Utility bath volume	280 l						
Dimensions: Length / Width / Height	2400 / 980 / 900 mm						
Main power supply	230 V / 50 Hz						

● – standard ○ – optional



ALFA 10

ALFA 20

ALFA 30
DeluxeALFA 40
Deluxe

ALFA 50

ALFA 60
DeluxeALFA 70
Deluxe



DEVELOPED
FOR FREQUENT
USE




HYDROTHERAPY

DELTA


FULL-IMMERSION BATHS WITH LOWER WATER VOLUME & OPERATIONAL COSTS



Hydromassage
jets



Air massage jets



Carbonated
water massage



Aromatherapy
and bath
additives



Colour therapy

FOR ALL TYPES OF MEDICAL FACILITIES & SPAS

ERGONOMICALLY-SHAPED

17



MAIN FEATURES & OPTIONS

- ◇ Utility bath volume of only 220 litres
- ◇ Easily-accessible even for people with lower mobility
- ◇ Comfortable operation using the colour touch display or buttons
- ◇ Ergonomic underwater massage
- ◇ Multi-zone therapy system
- ◇ 'Hot-keys' for frequently used therapies
- ◇ Pre-set therapy protocols
- ◇ Suitable for different kinds of water: Fresh, salt, carbonated or re-circulated
- ◇ Automatic disinfection programme
- ◇ High-quality acrylate bath surface for easy cleaning
- ◇ LED lighting of bath exterior
- ◇ Developed for reliability and long life

HYDROMASSAGE

- ◇ Up to 4 water jets and 38 micro-jets
- ◇ Up to 4 hydromassage zones:
 - Soles of the feet
 - Lower limbs, calves and knees
 - Buttocks, thighs and hips
 - Back, upper limbs and shoulders

AIR MASSAGE

- ◇ Up to 122 air massage jets
- ◇ Up to 2 air massage zones:
 - Lower limbs
 - Hips and buttocks

Model	DELTA 10	DELTA 20	DELTA 30	DELTA 40	DELTA 50 Deluxe
Specifications			DELTA 30 Deluxe	DELTA 40 Deluxe	DELTA 50 Deluxe
HYDROMASSAGE			●	●	●
Number of hydromassage jets					4
Number of hydromassage micro-jets			38	38	38
Hydromassage of the following body zones:			1 zone	4 zones	4 zones
Soles of the feet			1	1	1
Lower limbs				2	2
Outer calves and knees				3	3
Inner calves and knees				4	4
Buttocks					
Outer thighs and hips					
Lower and upper back					
Upper limbs and shoulders					
Regulation of hydromassage intensity			●	●	●
Hydromassage with added air			●	●	●
AIR MASSAGE		●		●	●
Number of air massage jets		122		122	122
Air massage of the following body zones:		1 zone		1 zone	2 zones
Lower limbs		1		1	1
Hips and buttocks					2
Regulation of air massage intensity		●		●	●
UNDERWATER MASSAGE	○	○	○	○	○
CARBONATED WATER BATH (CO₂ saturator)	●	○	○	○	○
STANDARD SPECIFICATIONS & OPTIONS					
Water level indicator with acoustic signaling	●	●	●	●	●
High water temperature indicator	●	●	●	●	●
Slide-in shower	●	●	●	●	●
Thermostatic mixing tap	○	○	○ (Deluxe)	○ (Deluxe)	●
Saltwater option / re-circulated water option	○	○	○	○	○
Foot support	○	○	○	○	○
External lighting	○	○	○	○	○
DELUXE LINE SPECIFICATIONS & OPTIONS					
LCD touch display 5.7"			● (Deluxe)	● (Deluxe)	●
Water filling & bath temperature indicator			● (Deluxe)	● (Deluxe)	●
Automatic bath filling and emptying			● (Deluxe)	● (Deluxe)	●
Thermostatic mixing tap			● (Deluxe)	● (Deluxe)	●
'Hot-keys' for frequently used programmes			● (Deluxe)	● (Deluxe)	●
Automatic drying of air massage jets				● (Deluxe)	●
Colour therapy option (Chromatherapy)			○ (Deluxe)	○ (Deluxe)	○
Music option (Built-in sound system)			○ (Deluxe)	○ (Deluxe)	○
Automatic disinfection programme			○ (Deluxe)	○ (Deluxe)	○
TECHNICAL SPECIFICATIONS					
Utility bath volume	220 l				
Dimensions: Length / Width / Height	2370 / 890 / 770 mm				
Main power supply	230 V / 50 Hz				

● – standard ○ – optional



DELTA 10

DELTA 20

DELTA 30
DeluxeDELTA 40
DeluxeDELTA 50
Deluxe

SPECIAL
ATTENTION GIVEN
TO EVERY SINGLE
COMPONENT



